



ALUMINA

30 mi.
Range

9.9 kg
Heaviest Part

21 st.
Weight Limit

The Alumina range offers the innovation of Lithium battery technology to a lightweight aluminium chassis, making it the easiest to lift and longest range travel scooter in our range.



Lithium-ion
Battery Technology



Quick release T-bar handle



Super soft tyres & full suspension



Front LED light



Rotating padded seat



Tiller charging socket



Dual controls & soft hand grips



Easily disassembles

Available in:



Charcoal



Grey

Motion

Specifications

	Air	Standard	Pro
Range	10 miles	15 miles	30 miles
Heaviest Part	9.9 kg	9.9 kg	9.9 kg
Top Speed	4 mph	4 mph	4 mph
Seat Width	45 cm	45 cm	45 cm
Wheels	19 cm	19 cm	19 cm
Battery Watt Hours	240 Wh	480 Wh	960 Wh
Total Weight	31.2 kg	33 kg	35 kg
Battery Weight	1.7 kg	3.5 kg	5.25 kg
Weight Limit	21 stone	21 stone	21 stone
Battery	10ah	20ah	40ah
Motor	270w	270w	270w
Ground Clearance	5 cm	5 cm	5 cm
Max. Gradient	6°	6°	6°
Turning Radius	125 cm	125 cm	125 cm

Accessories



Delta handlebars



Air safe battery



In-car charger

