

 **TOPRO**

Premium rollators



INT

made in Norway



We continue to expand our rollator and entire product portfolio

For several years now, TOPRO has continually updated, or added to, its rollator line-up – e. g. the Olympos ATR (All-Terrain-Rollator) or Pegasus, our lightweight carbon model. This time around, we are proud to present you with the TOPRO Hestia, our all-new indoor-rollator with never-before-seen features and benefits (for more details, turn to pages 6 to 10).

As the market leader for pioneering mobility solutions we don't stop there: our product portfolio now includes electric mobility scooters as well as AssiStep, a truly innovative, mechanical "rollator" for stairs (stair walker) – please refer to our website or the respective brochures for more information.

Developing and producing quality mobility aids has been our focus for over 30 years. We take pride in what we do and are happy to see how our products improve our customers' daily lives in the here and now.

Looking to the future, beyond mere product benefits, the environment also plays a big role in our thought process, with sustainability playing a central role: Our factory in Norway uses hydro-power and environmentally friendly aluminium to minimize the carbon footprint. Our injection-moulding plant uses clean materials that are recyclable to contribute to the circular economy. We have also implemented a "zero

waste"-principle to ensure no plastic contaminates the natural environment for the benefit of future generations.

Our goal, through our products, is to make your life easier. Please stay healthy, active – and „keep on moving“!



Asbjørn Eskild
CEO TOPRO Group

Advantages of the TOPRO rollator family

Whether it's product features or design – we have continually improved our rollators for over 30 years.

We are convinced our products provide you with an added value – by combining smart functionality and modern design with high quality and safety.

Unique memory function
Saves the preferred preset handle height, for easy re-use. (Standard on Troja 5G, Olympos ATR and Troja Original)

Easy wheel change
The interior braking system (TOPRO IBS) in combination with the TOPRO Quick Release system for easy wheel change, allows the use of different tyre treads. (Standard on TOPRO Pegasus, Troja 5G, Olympos ATR and Odysseé)



Ergonomic handles
The TOPRO Ergo Grip ensures correct walking posture and comfortable armrests while seated. Together with the soft tyres it provides double shock absorbing effect for stable driving comfort. (Standard on TOPRO Pegasus, Troja 5G, Olympos ATR and Odysseé)

Easy to fold
Simply pull the strap on the seat with only one hand. Practical to transport, stays firmly on its wheels and takes up little space when stowed.

Tilt function and edge guard
Enables the user to easily negotiate thresholds and kerbs, while keeping both hands on the rollator for safety.

Find the rollator that fits you best:
Our rollator guide assists you with selecting the right rollator for your individual needs:



TOPRO PEGASUS

CARBON FOR LIGHTWEIGHT, TOPRO FOR QUALITY



Modern, slender and lightweight carbon rollator with functionality and elegant design. This design combines to guarantee the highest level of quality, safety and stability at a weight of just 6.1 kg*.

Tried-and-tested features have been further developed:

- User-friendly method for setting the height of the handles using a **numeric scale**.
- **Partly integrated brake cables** and modern reflective panels on the frame.
- Elegant detachable shopping bag, which, with its retracted position provides **increased forward view to improve user safety**.
- Easy handling and manoeuvrability.
- **Reflectors** on the frame, lock wheels and shopping bag provide **360 degree visibility for safety in the dark**.
- **Shock-absorbing tyres** and **ball bearings** in wheels and forks for a comfortable ride.

Available in colour:

-  White
-  Steel blue
-  Wine red - NEW!

Measurements, colour options and other product details are available on page 18.



Learn more
TOPRO Pegasus

* Weight including shopping bag 6.6 kg.

“TOPRO Pegasus represents a new trend in the development of mobility aids where design, functionality and safety are combined to meet users’ expectations.”



made in Norway



HESTIA

FEEL RIGHT AT HOME



Stylish, easy-to-fold space-saving indoor rollator – offering you that extra bit of self-reliance and comfort at home. With a slew of clever features, the Hestia enables you to move about safely, get up and sit down more easily as well as to transport meals and other items securely.

You can look forward to these innovative features:

- Integrated, extendable **rise-up handles** provide reliable support for getting up and sitting down. The handles retract for unobstructed walking.
- **Multi-purpose tray** can easily be lifted and carried with one hand with its retractable handle, to safely transfer food and drinks.
- The **slidable tray** can be moved to the walker side, making it easier to reach when used as a table while seated.
- A **cup holder** and a **small groove** for your smartphone or tablet come standard. Clever: If you place your loudspeakers in the holder, their volume gets amplified!
- **Handrail** with a continuous brake lever allows for variable grip positions and safe handling while manoeuvring the compact Hestia.

Available in:



Silver

Measurements, colour options and other product details are available on page 19.



Learn more
TOPRO Hestia



“With Hestia, I can manage independently on my own at home and enjoy the beauty of everyday moments.”

| Karen



made in Norway

Hestia: The new benchmark for indoor rollators

The TOPRO Hestia is ready to go – no assembly required
It means you can get moving and use these practical details right away:



TOPRO Ergo Handle – Ergonomic and height adjustable, the one-piece handle-bar allows variable grip positions as well as single-handed operation.



Continuous Brake Lever – Safe and simple braking from different grip positions. Integrated parking brake.



Slidabe Tray – Tray slides back and forth for use as an easy-access table while seated. Integrated cup holder.

Carry with one Hand – Tray is easily detachable and balances itself when carried with one hand.



Basket – Discreet storage for personal items underneath tray.



Rise-up Handles – Retractable rise-up handles to support getting up from a seated position.



Tilt-function – To overcome small obstacles (e. g. thresholds)



TOPRO Hestia is both a sturdy and compact, lightweight indoor rollator, easy-to-maneuvre around your home. When you don't need it, it folds for compact storage (including the tray) and stands on its own, with a small footprint.

TROJA 5G

TO MAKE YOU FEEL BETTER



Handling comfort and stability – in true TOPRO-style. The TOPRO Troja 5G sets itself apart with its modern, slender design that goes hand-in-hand with a high level of quality and a secure driving experience.

It's the little things that make the difference – for improved usability and comfort:

- User friendly method for setting the height of the handles using a **numeric scale**.
- **Shock-absorbing tyres** provides a comfortable ride and smooth suspension.
- **Bell** to attract attention when necessary. Easy and safe to operate from the handle.
- Easy handling and manoeuvrability due to **ball bearings in forks and wheels**.

- **Detachable shopping bag** with reflective band, zipper and shoulder strap. Closeable and holds up to 10 kg. Two inner pockets for small items.
- Easy to fold, stays firmly on its wheels and takes up little space when stored.

Measurements, colour options and other product details are available on page 20.

Model M & S available in on-trend matt coated colours:



Learn more
TOPRO Troja 5G



“I have had a TOPRO rollator for a year now and it has transformed my life. I used to have a stick before, but once I got the rollator, I found I could walk so much better.“

| Audrey from the UK



OLYMPOS ATR

THE ALL-TERRAIN ROLLATOR



With large front wheels the Olympos ATR provides reliable and comfortable support on uneven roads and cobblestones as well as on sand or gravel foundations. The perfect choice for outdoor activities.

Innovative features to provide maximum support for outdoor use:

- **Wider design** to ensure stability in an outdoor environment.
- **Large and shock-absorbing tyres** to enable a comfortable ride and maintain stable and smooth manoeuvrability on rough terrain.
- **Detachable spacious water-resistant shopping bag** in on-trend grey fabric print. With reflective band and zipper. Closeable and holds up to 15 kg.

- User friendly method for setting the height of the handles with **numeric scale**.
- Easy to ride and manoeuvre due to **ball bearings in forks and wheels**.

Measurements, colour options and other product details are available on page 22.

Model M & S available in on-trend matt coated colours:



The ergonomic handle (TOPRO Ergo Grip) makes it easier and safer to get up from the sitting position. Together with the soft tyres it provides double shock absorbing effect for stable driving comfort.



Learn more
TOPRO Olympos ATR



made in Norway



TROJA ORIGINAL

THE EASY WAY TO GET AROUND





TOPRO's basic rollator with state-of-the-art design and high level of quality and features. The slender and modern design is infused with the functionality needed for a safe walk and an active everyday life.

The functionality is based on TOPROs high and continuous quality over many years:

- **Lightweight** at only 6.8 kg, and with a **weight bearing capacity** of 165 kg.
- Innovative **sideways folding technology** that is easy to operate with one hand by pulling the strap on the seat.
- **Ball bearings in forks and wheels** to ensure an easy ride and manoeuvrability.
- **Practical net basket** for easy access to items. **Reflective band** for safety and visibility in the dark.
- Available in three different sizes.
- Valued by users for being easy to handle and operate, having great manoeuvrability which makes it a trustworthy companion in everyday situations.
- Numeric scale for easy height adjustment of the handles.

Model M, S and X available in:

-  Silver
-  Wine red

Measurements, colour options and other product details are available on page 21.



Accessories



Crutch holder
can be fitted on either side.

Troja Original / 5G,
Olympos ATR,
Hestia – grey Art.No. 815358
TOPRO Pegasus Art.No. 815236
Olympos Art.No. 814040
Odysseé Art.No. 814457

*TOPRO Pegasus, Troja 5G,
Olympos ATR, Troja Original,
Odysseé, Troja Neuro, Hestia*



Tray
Snaps in securely, easy to remove
and clean, carries up to 5 kg

TOPRO Pegasus,
Troja Original / 5G,
Odysseé Art.No. 814728
Olympos ATR M Art.No. 815212

*TOPRO Pegasus, Troja 5G,
Olympos ATR, Troja Original,
Odysseé*



Back support with padding
ergonomically height adjustable
according to preferred seating
position. Available in two
lengths.

Short 67 cm Art.No. 815270
Long 75 cm Art.No. 815272

*TOPRO Pegasus, Troja 5G, Olympos
ATR, Odysseé, Troja Neuro*



**LED lamp with clip on and USB
plug for safety in the dark.**
Chargeable.

Art.No. 815260

*TOPRO Pegasus, Troja 5G,
Olympos ATR, Troja Original,
Odysseé, Troja Neuro*



Back support with padding
ergonomically height adjustable
according to preferred seating
position. Available in two
lengths.

Short 68 cm Art.No. 815271
Long 76 cm Art.No. 815273

Troja Original



Umbrella with attachment
keep dry under the umbrella in
rainy weather or use it as a parasol.

Black Art.No. 814068
Red Art.No. 814099

*TOPRO Pegasus, Troja 5G,
Olympos ATR, Troja Original,
Odysseé, Troja Neuro*



Additional accessories can be found on our website.

Accessories



made in Norway



Technical data for TOPRO Pegasus



Recommended user height	150–200 cm
Maximum weight of user	130 kg
Maximum capacity shopping bag	10 kg / 15 l
Handle height	78–97.5 cm
Handles	Ergonomic, several holding positions (TOPRO ErgoGrip)
Grip distance driving brake	7 cm
Width between center points of handles	47 cm
Seating width	43.5 cm
Seat height	63 cm
Turning circle	84 cm
Maximum height	101 cm
Maximum width	64 cm
Maximum length	65 cm
Diameter wheels	20 cm
Width wheels	3.6 cm
Length when folded	65 cm
Width when folded	27 cm
Height when folded	80.5 cm
Dimensions folded (l × w × h)	65 × 27.3 × 80.3 cm
Product weight (excl. bag)	6.1 kg
Product weight (incl. bag)	6.6 kg
Area of use	Indoors & outdoors
Article No.	815020 + item No.

Available in colour:



Technical data for Hestia



Recommended user height	156–196 cm
Maximum weight of user	110 kg
Maximum load storage net	5 kg
Handle height	81–100 cm
Handlebar	Ergonomic, several holding positions (TOPRO ErgoBar)
Grip distance driving brake	7.5 cm
Width between center points of handles	37.5 cm
Maximum load tray	5 kg
Height of tray position	68.3 cm
Tray size	35 × 29 cm
Turning circle	78 cm
Maximum height	100 cm
Maximum width	55.5 cm
Maximum length	63 cm
Diameter wheels	11.2 cm
Width wheels	2.4 cm
Length when folded	25 cm
Width when folded	55.5 cm
Height when folded	88.5 cm
Dimensions folded (l × w × h)	25 × 55.5 × 88.5 cm
Product weight	6.9 kg
Area of use	Indoors
Article No.	815500-119010

Available in colour:



Technical data for Troja 5G



Size	M (medium)	S (small)
Recommended user height	150–200 cm	135–170 cm
Maximum weight of user	150 kg	125 kg
Maximum capacity shopping bag	10 kg / 15 l	10 kg / 15 l
Handle height	78–100 cm	69–86 cm
Handles	Ergonomic, several holding positions (TOPRO ErgoGrip)	
Grip distance driving brake	7 cm	7 cm
Width between center points of handles	47 cm	47 cm
Seating width	45.5 cm	45.5 cm
Seat height	63 cm	54 cm
Turning circle	84 cm	84 cm
Maximum height	103.5 cm	89.5 cm
Maximum width	60 cm	60 cm
Maximum length	65 cm	65 cm
Diameter wheels	20 cm	20 cm
Width wheels	3.6 cm	3.6 cm
Length when folded	65 cm	65 cm
Width when folded	24 cm	24 cm
Height when folded	79.5 cm	73 cm
Transport dimensions (l × w × h) (excl. wheels and shopping bag)	65 × 24 × 73 cm	65 × 24 × 66 cm
Transport weight (excl. wheels and shopping bag)	5.6 kg	5.4 kg
Product weight (incl. wheels and shopping bag)	7.8 kg	7.6 kg
Area of use	Indoors & outdoors	
Article No.	814600 + item No.	814601 + item No.

Model M & S are available in on-trend matt coated colours:

-  Silver
-  Sand
-  Yellow
-  Wine red
-  Purple
-  Blue
-  Turquoise
-  Black

Technical data for Troja Original



Size	M (medium)	S (small)	X (extra small)
Recommended user height	150–200 cm	130–170 cm	120–160 cm
Maximum weight of user	165 kg	165 kg	125 kg
Maximum capacity basket	10 kg / 20 l	10 kg / 20 l	5 kg / 11 l
Handle height	78–100 cm	69–86 cm	64–79 cm
Handles	Anatomical		
Grip distance driving brake	6.5 cm	6.5 cm	6.5 cm
Width between center points of handles	47 cm	47 cm	37 cm
Seating width	45.5 cm	45.5 cm	36 cm
Seat height	63 cm	54 cm	51 cm
Turning circle	84 cm	84 cm	79 cm
Maximum height	101 cm	87 cm	80 cm
Maximum width	60 cm	60 cm	51 cm
Maximum length	65 cm	65 cm	65 cm
Diameter wheels	20 cm	20 cm	20 cm
Width wheels	3.6 cm	3.6 cm	3.6 cm
Length when folded	65 cm	65 cm	65 cm
Width when folded	24 cm	24 cm	24 cm
Height when folded	80 cm	74 cm	66 cm
Dimensions folded (l × w × h)	65 × 24 × 80 cm	65 × 24 × 74 cm	65 × 24 × 65 cm
Product weight (incl. wheels and basket)	6.8 kg	6.6 kg	6.1 kg
Area of use	Indoors & outdoors		
Article No.	815400-110000	815410-111000	815420-112000

Model M, S and X are available in on-trend matt coated colour:

-  Silver


Technical data for Olymos ATR



Size	M (medium)	S (small)
Recommended user height	150–195 cm	140–170 cm
Maximum weight of user	150 kg	150 kg
Maximum capacity shopping bag	15 kg / 12 l	15 kg / 10 l
Handle height	76–95.5 cm	69–86.5 cm
Handles	Ergonomic, several holding positions (TOPRO ErgoGrip)	
Grip distance driving brake	7 cm	7 cm
Width between center points of handles	49 cm	50 cm
Seating width	47 cm	47 cm
Seat height	60 cm	54 cm
Turning circle	97 cm	97 cm
Maximum height	99 cm	90 cm
Maximum length	77 cm	77 cm
Maximum width	68 cm	68 cm
Diameter wheels	28 cm / 23 cm	28 cm / 23 cm
Width wheels	4.2 cm	4.2 cm
Length when folded	77 cm	77 cm
Width when folded	27 cm	27 cm
Height when folded	87.5 cm	84 cm
Dimensions folded (l x w x h)	77 x 27 x 87.5 cm	77 x 27 x 85 cm
Product weight (incl. wheels and shopping bag)	8.9 kg	8.7 kg
Area of use	Outdoors	
Article No.	814300 + item No.	814305 + item No.

Model M & S are available in on-trend matt coated colours:

-  Silver
-  Sand
-  Yellow
-  Wine red
-  Purple
-  Blue
-  Turquoise
-  Black



The future is now

Choose a TOPRO rollator and help ensure a sustainable future.
Read more about our contribution on topromobility.com

Global responsibility

Sustainable materials and manufacturing

Clean energy

Sustainability in focus

We embrace the future by reducing our environmental footprint and contributing more to the circular economy, with implementing new accessory packaging made from materials that are easy to recycle.



Premium rollators from Norway

Manufacturer:



TOPRO Industri AS

Rambekkvegen 5

2816 Gjøvik

Norway

Tel. (+47) 61 13 46 00

info@topromobility.com

www.topromobility.com



Learn more



topromobility.com

Designed and manufactured in Norway.

TOPRO vouches for top quality
and offers a 7-year warranty.

Setup and variations in fixed accessories for
the rollators may differ between countries.

Quality Assurance Standards

EN ISO 9001 – quality

EN ISO 14001 – environment

EN ISO 11199-2 – rollator

EN ISO 10993-5 – biocompatibility

made in Norway



Dealer